


Manufacturer	Sartorius
Model number	Sterisart Universal Gen 4
Serial number	740200022
Package version	09-09-01.00.16
Task name	Pomper, scanner et enregistrer
Device ID 1	SALF SPA
Device ID 2	LAB MICROBIOLOGIA

User	Administrateur
Workflow configuration	Pump, Scan & Record
Overpressure response	Manual regulation
Start of process	2025-06-25 14:43:53
Sample ID	809008035



Canister Type	16466-----GS
---------------	--------------

Run 1:	Start: Speed: 50. Duration: Off. Pump started: 2025-06-25 14:44:09. User stopped the pump at 0:12 min:sec. Pump stopped: 2025-06-25 14:44:21. Result: Average speed: 50. Duration: 0:11 min:sec.
End of process	2025-06-25 14:44:21

Date/time of report end	2025-06-25	14:44:44
Signat.		